
USB cable connection driver installation SOP

- In general, when pettrust is connected to a computer via a USB cable, the computer will download and install the driver itself.
- If the client's computer cannot be transferred or the COM cannot be seen by the device administrator, the following steps can be performed.
- The following steps are based on Windows 7

1. Right click on "My Computer" and click on the content.

2. Confirm the computer system version (such as Windows 7 on the left) and the system type (32 bits on the left).

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3. Select the corresponding driver installation, confirm the Windows 7_32 bit as above, and select "CP210x_Driver_Win7.8.10" as shown on the left.

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
4. Select the application installation of "CP210xVCPInstaller_x86 (32 bit)" as shown on the left.

5. Follow the instructions for driver installation.

6. After the installation is complete, connect to pettrust, enter the device administrator, and click “COM&LPT” to confirm whether there is “Silicon Labs CP210x USB to UART Bridge (COM4)” display. If yes, the connection is successful.